

# ProClear™

## ORTHODONTIC ACCELERATOR INSTRUCTIONS FOR USE



ProClear™ Orthodontic Accelerator provides low frequency vibrations (30Hz) designed to help alleviate the pain and discomfort associated with orthodontic braces and invisible aligners, while promoting faster results.

### **KIT CONTENTS:**

- Orthodontic Accelerator & Mouthpiece
- USB Cable
- Charger & Wall Socket Plug
- Charger Cradle

## CHARGING INSTRUCTIONS:

- Place orthodontic accelerator in charger cradle
- Insert USB cable into charger cradle
- Plug into power source
- Charge device for 2 hours
- The button will flash blue when a recharge is needed

## OPERATION:

- Attach the included mouthpiece to orthodontic accelerator
- Press the device button to power on
- Place the device in your mouth and begin treatment using the recommended times below:
  - 10 minutes in the morning
  - 10 minutes in the evening

## HOW TO CLEAN:

- Remove the mouthpiece from the device
- Clean mouthpiece with water

**FAST RESULTS:** Using the orthodontic accelerator for just 20 minutes a day helps accelerate the movement of teeth and reduces the overall treatment time.

Visit us online for additional help, including tips and tutorials please visit

[www.ProClearAligners.com/Support](http://www.ProClearAligners.com/Support)